



## HUMAN PLASMA, RHEUMATOID FACTOR (RF) POSITIVE, TITER >2,500 IU/ML

Catalog Number: 0101000550  
Source: Human  
Storage: -20° C  
Expiration: 3 years from DOM

Lot Number: Lot Specific  
Fill: Single Donor(s)  
Preservative: None  
Filtration: None

Specification Sheet

### Product Specifications

Analyte	Result	Unit of Measure
Rheumatoid Factor (RF)	>2,500	IU/ml

### Testing

This material has been tested for HIV 1/2, HIV-1 (NAT), HBsAg, HBV (NAT), HCV, HCV (NAT), and Syphilis by FDA approved methods and found to be NON-REACTIVE.

### Caution

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulations.

### Warranty

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.