

BETA-2 MICROGLOBULIN

Catalogue Number: A1010

Lot Number: *Lot Specific*

Source: Human Urine

Perservative: None

Storage: -20° C

Product Performance Specifications

Test	Result
Appearance	White, Lyophilized Powder
Form	Supplied lyophilized from 0.02 M NH ₄ HCO ₃ (may contain traces of buffer salts)

NOTE:

Immunological identity and β -2-M concentration determined by Radical Immunodiffusion (RID). Electrophoretic and Chromatographic data of the major component is consistent with molecular weight of β -2-M (approximately 11,800).

TESTING:

This material has been tested for HIV 1/2, HIV-1 (NAT), HBsAg, HCV, and HCV (NAT) by FDA approved methods and found to be NON-REACTIVE.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

