

HUMAN CHORIONIC GONADOTROPIN (HCG)

CATALOG NUMBER: A1050 LOT NUMBER: TBD

Storage:	2 - 8°C	Expiration Date:	<Lot Specific>
Source:	Human	Manufacture Date:	<Lot Specific>
Appearance:	Lyophilized Powder	Clarity of solution:	Clear
Color of solution:	Colorless		

Product Performance Specifications

Test	Specifications	Unit of Measurement
Solubility (5mg/ml)	≤ 2 minutes	N/A
Identification HCG	Positive	N/A
Water	≤5	%
PH	6.0 – 8.0	N/A
Potency HCG	≥ 1,500	IU/mg
Sterility	Sterile	N/A
Acute Toxicity (USP)	No symptoms of toxic reaction	N/A
Estrogenic activity (USP)	Complies to USP	N/A
Residual solvents	Complies to USP	N/A
Bacterial Endotoxins (USP)	≤ 0.03	USP EU / I.U. HCG

TESTING:

This material has been tested for HIV 1/2, HBsAg, and HCV by FDA approved methods and found to be NON-REACTIVE.

CAUTION:

In-Vitro use only. Products are intended for research or manufacturing of non-injectable products only. This material was tested as stated above in the testing section. Nevertheless, products of human origin should be considered potentially infectious and handled according to universal precautions. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury or illness resulting from the use and/or handling of the product. The customer assumes responsibility for the proper disposal of all materials in accordance with Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

