



27625 Commerce Center Dr  
Temecula, CA 92590

---

## ***REPRESENTATIVE CERTIFICATE of ANALYSIS***

---

PRODUCT:	HUMAN C-REACTIVE PROTEIN (CRP)
CATALOG #:	A1084
LOT #:	<i>Lot Specific</i>
SOURCE:	Human Fluids
FORM:	Liquid
BUFFER:	20 mM Tris, 280 mM Sodium Chloride, 5 mM Calcium Chloride
PRESERVATIVE:	0.09% Sodium Azide
pH:	<i>Lot Specific</i>
CRP CONTENT:	<i>Lot Specific</i>
TOTAL PROTEIN:	<i>Lot Specific</i>
PURITY:	<i>Lot Specific</i>
STORAGE:	2-8° C (REFRIGERATE, DO NOT FREEZE)
EXPIRATION:	<i>Lot Specific</i>

**TESTING:** This material has been tested for HIV 1/2, HIV-1 (NAT), HBsAg, HBV (NAT), HCV and HCV (NAT) by FDA approved methods and found to be NON-REACTIVE.

**CAUTION:** Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulations.

**WARRANTY:** Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.