

HEPATITIS B SURFACE ANTIGEN SUBTYPE AD

Catalogue Number: A1280-AD

Lot Number: *Lot Specific*

Source: Human

Storage: 2-8° C (REFRIGERATE, DO NOT FREEZE) Fill: Single Donor

Product Performance Specifications

Test	Result
Form	Liquid (0.2 micron filtered)
Buffer	20 mM Phosphate buffer at pH 7.4 containing 150 mM NaCl 5% sucrose and 15 mM NaN ₃

TESTING:

This material has been tested for HIV 1/2, HIV-1 (NAT), HBsAg, HBV (NAT), HCV, HCV (NAT), and RPR by FDA approved methods and found to be NON-REACTIVE.

CAUTION:

NOT FOR USE IN HUMANS OR CLINICAL DIAGNOSIS. These products are intended for research or manufacturing of NON-INJECTABLE products only.

"PRODUCT TESTED POSITIVE FOR HEPATITIS B and the risk of transmitting HEPATITIS B is present. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation."

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

