

CANINE SERUM

Catalogue Number: AS1020-A

Lot Number: Lot Specific

Source: Normal Non-Specific Canines

Perservative: None

Storage: -20° C or Colder

Product Performance Specifications

| | Bleed Number | Bleed Date | Age | Sex | Race | Volume (mL) |
|---|--------------------|------------|-----|-----|------|-------------|
| 1 | Donation Specific* | * | * | * | * | * |

NOTE:

Although this material was obtained from canines that were determined to be healthy by a doctor of veterinarian medicine (D.V.M.) for an extended amount of time prior to collection, please handle these biologicals as if they are capable of transmitting disease. Serum was obtained aseptically from spontaneously clotted whole blood from healthy, adult domestic canines using venipuncture technique. All canines used are currently conditioned for pharmaceutical research laboratory use and were fasted for a minimum of 8 hours prior to collection.

TESTING:

This material has NOT BEEN VIRAL TESTED.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

