

BOVINE CHOLESTEROL CONCENTRATE

Catalogue Number: BF3090

Lot Number: *Lot Specific*

Source: Bovine Serum
Storage: -20° C or Colder

Perservative: None

Product Performance Specifications

Test	Result
Appearance	Light yellow to brown liquid

NOTE:

The animal source of this product originates in the United States of America and is collected at a USDA licensed establishment. The animals receive ante- and post-mortem health inspections at the abattoir by a US FSIS inspector. The animals are apparently healthy and free from infectious and contagious diseases and injurious parasites.

TESTING:

This material has NOT BEEN VIRAL TESTED.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

