

HUMAN PLASMA, LI HEPARIN, SARS

Catalogue Number: D2008

Lot Number: Lot Specific

Source: Human

Perservative: None

Storage: -20° C

Fill: Single Donors

Product Performance Specifications

	Bleed Number	Bleed Date	Age	Sex	Race	Volume (mL)
1	Donation Specific*	*	*	*	*	*

NOTE:

This material was collected under an IRB approved protocol.

This material tested positive for IgG antibodies.

This material has been tested using FDA approved assays. Only the analyte(s) described above have been tested and passed internal QC specifications with the designated assay method. Values obtained from other manufacturers' assay methods may not result in corresponding values. This difference does not invalidate results supplied.

Use only for the manufacture of unlicensed, non-injectable diagnostic products not subject to license by the FDA under Section 351 of the Public Health Service Act.

The supplier does not determine expiration dates on bulk human sera for use in diagnostic manufacturing products, not knowing the final output of how the material is used by the end user. The expiration date must be determined by the end user.

TESTING:

This material has NOT BEEN VIRAL TESTED.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation. This material should be handled at Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimens in the CDC/NIH manual "Biosafety in Microbiological and Biomedical Laboratories.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

