

HUMAN PLASMA, HBSAG AD POSITIVE

Catalogue Number: D4141

Lot Number: *Lot Specific*

Source: Human

Perservative: None

Storage: -20° C

Fill: Single Donor

Product Performance Specifications

	Bleed Number	Bleed Date	Age	Sex	Race	Volume (mL)
1	Donation Specific*	*	*	*	*	*

NOTE:

Material was tested for presence of Ad subtype using HBsAg Subtype EIA Kit. The donor of each unit was tested at the donor level, not for every donation.

TESTING:

This material has been tested for HIV 1/2, HCV, HIV-1 (PCR) and HCV (PCR) by FDA approved methods and found to be NON-REACTIVE. This material has also been tested for Syphilis per FDA guidelines and found NEGATIVE. PRODUCT TESTED POSITIVE FOR HEPATITIS B and the risk of transmitting HEPATITIS B is present.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

