

GOAT ANTI-HUMAN IGG FC, NEPHELOMETRIC GRADE ANTISERA

Catalogue Number: NA2000

Lot Number: *Lot Specific*

Source: Goat

Perservative: 15 mM NaN3

Storage: 2-8° C (REFRIGERATE, DO NOT FREEZE)

Product Performance Specifications

Test	Result
Appearance	Liquid
Form	Liquid (Affinity purified Polyclonal Antibody, delipidated, and 0.2 micron filtered)
Filtration	0.2 micron filtered
Buffer	20 mM Tris containing 250 mM NaCl, pH 7.4
Specificity	Immunoelectrophoresis gives single precipitin arc specific for human IgG Fc, with normal human serum

TESTING:

This material has NOT BEEN VIRAL TESTED.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

