

## GOAT ANTI-HUMAN TROPONIN I, ANTIBODY

**Catalogue Number:** NA2023

**Lot Number:** Lot Specific

Source: Goat

Perservative: 0.09% NaN3

Storage: -20° C

United States

### Product Performance Specifications

Test	Result
Form	Polyclonal, Delipidated, Affinity Purified
Filtration	0.2 µm
Buffer	0.05 M Tris, 0.5 M NaCl, 0.01 M EDTA, pH 7.5
Specificity	By IEP, no reaction when tested against 2xPHS. Positive reaction against Troponin I.
Antigen	Purified Human Troponin I

#### NOTE:

If particulates are observed, filtration prior to use is recommended. The supplier does not determine expiration dates on this material for use in diagnostic manufacturing products, not knowing the final output of how the material is used by the end user. The expiration date must be determined by the end user.

#### TESTING:

This material has NOT BEEN VIRAL TESTED.

#### CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

#### WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

