

## Normal Human Urine, Male

Catalog Number: OH2000

Lot Number: Lot Specific

Storage: -20 °C  
Manufacture Date: Lot Specific  
Country of Origin: USA

Source: Male Human Urine  
Preservative: None  
Expiration Date: Lot Specific

### Product Performance Specifications

Test	Specification	Test	Specification
Specific Gravity	≥ 1.007 g/mL	Urine Nitrites	Negative
HCG	<5 mIU/mL	Urobilinogen	0.2-1.05 EU/dL
Urine Protein	Negative	Leukoesterase	Negative
Urine Glucose	Negative	Color	Yellow
Urine Ketones	Negative	WBC	0-4 HPF
Urine Bilirubin	Negative	RBC	0-2 HPF
Urine Blood	Negative	Urine Potassium	≤ 68 mmol/L
Odor	Urine	Urine Sodium	< 400 mmol/L
pH	5.5 – 9.0		

### Testing

Donors have been tested and found negative or non-reactive for HIV 1/2, HBsAg, HCV, HIV-1 (NAT), HCV (NAT) and RPR by FDA Approved Methods.

### Caution

In-vitro use only. Products are intended for research or manufacturing of non-injectable products only. This material was tested as stated above in the testing section. Nevertheless, products of human origin should be considered potentially infectious and handled according to universal precautions. The safe and responsible handling of the product listed above is the responsibility of the customer. We accept no liability for injury or illness resulting from the use and/or handling of the product. The customer assumes responsibility for the proper disposal of all materials in accordance with Federal and State regulations.

### Warranty

Golden West Diagnostics, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.



Golden West Diagnostics, LLC · 42191 Sarah Way, Temecula, CA 92590 · +1-951-699-0957