

## GOAT ANTI-MOUSE IGG, AFFINITY PURIFIED

**Catalogue Number:** SA1060(AP)

**Lot Number:** Lot Specific

---

**Source:** Goat

**Storage:** 2-8° C (REFRIGERATE, DO NOT FREEZE)

### Product Performance Specifications

Test	Result
Form	Nonsterile liquid. Phosphate buffered saline with 0.05% Sodium Azide. Filtered through a 0.2-micron membrane filter.
Filtration	0.2-micron
Specificity	Reacts with Mouse IgG gamma chains, as well as the light chains from all Mouse Immunoglobulin classes. No reactivity is observed to non-immunoglobulin Mouse serum proteins.
Immunogen	Purified Mouse IgG, whole molecule

#### NOTE:

Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

#### TESTING:

This material has NOT BEEN VIRAL TESTED.

#### CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

#### WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

