

GOAT ANTI-MOUSE IGG

Catalogue Number: SA1060

Lot Number: Lot Specific

Source: Goat

Perservative: None

Storage: -20° C

United States

Product Performance Specifications

Test	Result
Form	Nonsterile liquid. Filtered through a 0.45-micron membrane filter.
Filtration	0.45-micron
Specificity	Reacts with Mouse IgG gamma chain, as well as the light chains from all Mouse Immunoglobulin classes. Cross reactivity with immunoglobulins from other species is also observed.
Immunogen	Highly purified Mouse IgG isolated from pooled normal Mouse serum.

NOTE:

The supplier does not determine expiration dates not knowing the final output of how the material is used by the end user. The expiration date must be determined by the end user.

TESTING:

This material has NOT BEEN VIRAL TESTED.

CAUTION:

Products are intended for research or manufacturing of non-injectable products only. The safe and responsible handling of the products listed above is the responsibility of the customer. We accept no liability for injury resulting from their use or handling. The customer assumes responsibility for the proper disposal of all materials in accordance with proper Federal and State regulation.

WARRANTY:

Golden West BioSolutions, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material, and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.

